•  
Sample collection  
and submission guidelines  
•


*Paper Submission only for remaining kits distributed before April 2024.*

# Content of the kits received:


## 1. Collection form

## 2. Test request form

## 3. Physician's statement

# Content of the kits received:

**4.** ICF  
(Consent term)  
– 2 copies

**sanofi**

INFORMED CONSENT FORM TO TAKE PART IN THE SANOFI GLOBAL DIAGNOSTIC SUPPORT PROGRAM

Identification Data of the participant

Physician in charge: \_\_\_\_\_ ID No.: \_\_\_\_\_

Patient name: \_\_\_\_\_

Date of Birth: \_\_\_\_/\_\_\_\_/\_\_\_\_ Patient ID: \_\_\_\_\_

Legal guardian (if applicable): \_\_\_\_\_

You have been referred by your physician to take part in the Sanofi Global Diagnostic Support Program ("Program"), which is only liable for the allowance to conduct diagnostic exams requested by your physician and available in the Program Rules, subject to amendments with no prior notice. Any other required conduct, including, but not limited to, supplementary exams which are not available in the Program, or suspended temporarily, medical or other professional care, other medical services and provision of treatment are not Sanofi's liability, taking into account the Program's features.

This form must be read carefully and, should you have any questions, you can contact the physician responsible for your monitoring. After the following information has been clarified, we ask you to sign at the bottom of this document if you agree to use the services of the Program according to the terms below. There are two (2) copies. One is yours and the other is for the Laboratory hired by Sanofi (that is, Instituto Hermes Pardun SA in Brazil).

Consent to the terms below is a necessary condition for using the Program services. If you do not agree with the provisions of these terms, refusing to sign the document will make it impossible for the Laboratory hired by Sanofi to provide such services.

**I declare that I agree to enjoy the benefits of the Program freely and voluntarily and that I have been informed about the following topics:**

1. The exams are intended to assist in the research of the supposed diagnosis of pre-existing rare diseases set in the Sanofi Global Diagnostic Support Program, as requested by your physician, and **your participation in this Program is not conditioned to the payment of any amount, that is, you will not receive or pay any sums to participate in the Program.**
2. The objective of the Program is to perform laboratory analyses (enzymatic, molecular, genetic puncts, or others) using the appropriate and available techniques for diagnosis, through the collection of blood, urine, saliva, or other biological material, for the purpose of investigating the diagnostic suspicion free of charge, according to medical request.
3. The services provided under the Program are funded by Sanofi and comply with national and international quality guidelines, being considered the most suitable for the identification of pre-existing diseases.

global-icse-ssp-part-ing-v1-2023

**5.** Plastic envelope  
for return

**6.** Silica gel bag

**7.** External  
envelope

*Step 1:*

**Fill all the blanks in the forms within the kit.**


# Rare Disease Gene Panel Program

## Form 3: *Physician's Statement*

Fill with Physician's Information

Physician's signature

**sanofi**

**PHYSICIAN STATEMENT  
DIAGNOSTIC TESTING SERVICES BY  
SANOFI GLOBAL DIAGNOSTIC SUPPORT PROGRAM**

I, \_\_\_\_\_, e-mail \_\_\_\_\_,  
practicing at \_\_\_\_\_, department \_\_\_\_\_, city \_\_\_\_\_,  
country \_\_\_\_\_ hereby certify that I requested Sanofi to provide access to diagnostic testing without  
charge for the diagnosis of rare diseases.

The requested diagnostic testing will be performed by specialized laboratories retained by Sanofi under a service agreement. Sanofi is only coordinating and supporting this diagnostic service and is not otherwise involved in the diagnostic testing.

In order to honor this request, I, the undersigned, agree to adhere to the diagnostic testing services process (forms and systems) as outlined below:

1. use of Dried Blood Spot (DBS) cards + DBS sampling guidelines provided by laboratory contracted by Sanofi
2. use of a laboratory internal system to register the DBS sample, track its status (if applicable) and receive notifications when the analysis report is ready.
3. use of Informed Consent Forms (ICF), with the ICF template provided by laboratory contracted by Sanofi.

By signing this statement, I, the undersigned, confirm the following:

- the diagnostic testing is performed at my own request and any therapeutic decisions resulting from this testing are at my own responsibility.
- I will not seek nor accept reimbursement from a 3rd party for the services provided by Sanofi.
- I will not be paid or rewarded in any way by Sanofi to test patients and/or use the diagnostic services.
- I declare there is no conflict of interest between me (or any of my family members) and Sanofi.
- I guarantee that I will timely inform my institution of the requested services and the content of this statement.

\_\_\_\_\_  
Signature Date: \_\_\_\_/\_\_\_\_/\_\_\_\_

\_\_\_\_\_  
Professional title and full name

Your personal data will be processed in the context of laboratory testing. If you have questions on this processing, wish to receive a full privacy notice or exercise your access right or other statutory rights, please contact [maledjurni@slc.com.br](mailto:maledjurni@slc.com.br).


Fill with Physician's Name and  
Professional Title

# Rare Disease Gene Panel Program

## Form 4: Informed Consent Form (ICF)

Fill with Patient Information

**sanofi**

**INFORMED CONSENT FORM TO TAKE PART IN THE SANOFI  
GLOBAL DIAGNOSTIC SUPPORT PROGRAM**

**Identification Data of the participant** APPLY TAG B

Physician in charge: \_\_\_\_\_ ID No.: \_\_\_\_\_

Patient name: \_\_\_\_\_

Date of Birth: \_\_\_\_/\_\_\_\_/\_\_\_\_ Patient ID: \_\_\_\_\_

Legal guardian (if applicable): \_\_\_\_\_

You have been referred by your physician to take part in the Sanofi Global Diagnostic Support Program ("Program"), which is only liable for the allowance to conduct diagnostic exams requested by your physician and available in the Program Rules, subject to amendments with no prior notice. Any other required conduct, including, but not limited to, supplementary exams which are not available in the Program, or suspended temporarily, medical or other professional care, other medical services and provision of treatment are not Sanofi's liability, taking into account the Program's features.

This form must be read carefully and, should you have any questions, you can contact the physician responsible for your monitoring. After the following information has been clarified, we ask you to sign at the bottom of this document if you agree to use the services of the Program according to the terms below. There are two (2) copies. One is yours and the other is for the Laboratory hired by Sanofi (that is, Instituto Hermes Pardini SA in Brazil).

Consent to the terms below is a necessary condition for using the Program services. If you do not agree with the provisions of these terms, refusing to sign the document will make it impossible for the Laboratory hired by Sanofi to provide such services.

**I declare that I agree to enjoy the benefits of the Program freely and voluntarily and that I have been informed about the following topics:**

1. The exams are intended to assist in the research of the supposed diagnosis of pre-existing rare diseases set in the Sanofi Global Diagnostic Support Program, as requested by your physician, and **your participation in this Program is not conditioned to the payment of any amount, that is, you will not receive or pay any sums to participate in the Program.**
2. The objective of the Program is to perform laboratory analyses (enzymatic, molecular, genetic panels, or others) using the appropriate and available techniques for diagnosis, through the collection of blood, urine, saliva, or other biological material, for the purpose of investigating the diagnostic suspicion free of charge, according to medical request.
3. The services provided under the Program are funded by Sanofi and comply with national and international quality guidelines, being considered the most suitable for the identification of pre-existing diseases.

global-tcle-esp-port-ing-v1-2023

← Apply Tag A

- ⚠ Patient or guardian signature on last page. **(mandatory)**
- ⚠ 2 copies (one copy must be attached to the test request form and the other to be given to the patient).

*Step 2:*  
**Sample Collection Process.**


## Rare Disease Gene Panel Program

### Step 2:

### *Tips for Successful Sample Collection*

- Sufficient and homogeneous concentration of sample.
- Blood should be viewed in both sides of the filter.
- Sample could be drawn with needle or by capillary puncture.
- Important to fill all the paper filter circle with the blood.
- If the sample doesn't meet the quality standard, a new sample collection will be asked.

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

E-mail / Correo Electrónico

**Important!**  
The sample must be left drying for at least 3 hours.

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• Perform the collection by following the instructions for each sample on the right.

• Check the type of sample collected.

• Paste the tag A on the filter paper end the tag B on the test request form.

• Let it drying naturally, in a clean, horizontal and polished area, for at least 3 hours.

• Never use ways to accelerate drying (stove, sun, hair dryer, microwave, fan etc.).

• After the drying is completed, place the sample into the envelope form, seal and insert into protection plastic envelope.

• Keep it in room temperature (18°C to 25°C) until shipping to DLE Laboratory.

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

dle.com.br

*Step 3:*  
**Sample Packaging and Submission.**

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TRIAGEM FAMILIAR/FAMILY SCREENING/ESTUDIO FAMILIAR

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Place the dried blood spot card into the form envelope and seal it.

Insert the form envelope inside the plastic bag with silica gel and seal it.

The logo features a stylized DNA double helix on the left, composed of purple dots and lines. To its right, the text "RD Gene PP" is displayed in a bold, sans-serif font. "RD" and "PP" are in black, while "Gene" is in purple.

# RD Gene PP

Rare Disease Gene Panel Program

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**Proibida a reprodução total ou parcial deste material. Todos os direitos reservados ao Grupo Parfith. Total or partial reproduction of this material is forbidden. All Rights reserved to Grupo Parfith. Prohibida la reproducción total o parcial de este material. Todos los derechos reservados a Grupo Parfith.**